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JPY210,000 (JPY200,000 Tax excluded)

SAKAMOTO MODEL
CORPORATION

Basis of respiratory management Endotracheal intubation can be learned.

This model is a simulator with a lifelike texture, and various suction methods can be learned.

The suction of secretions not only in the oral/nasal cavities but also in tubes inserted via the oral/nasal cavities is possible.

Since the trachea has already been incised, suction with tracheal cannula insertion can also be performed.

In addition, this model incorporates bilateral lungs, allowing confirmation of the state of air delivery using an aspirator or a bag-valve-mask system.

Various types of suction training are possible

Suction of secretions in the cavity

Suction of secretions in the nasal cavity

Suction of secretions in inserted tubes


Suction of secretions in the tracheal cannula

SARACOTO MODEL
JAPAN

The head can be moved to the left and right.

Precautions

- ① Since the head may be deformed on contact with resins such as polyvinyl chloride, avoid its direct contact with resin-based products.
- ② The use of rubber gloves is recommended since the head easily becomes dirty due to the properties of its composing material. Once made, marks from ball-point or marker pens cannot be removed.
- ③ Slight stains on the head can be washed off using neutral detergents. (Do not use thinner because resin causes deformation.) Subsequently, wipe off water well and apply baby powder.
- ④ This product is a simulator allowing advanced training but does not guarantee similar effects in actual practice.
- ⑤ For improvement, the specifications of this product are subject to change without prior notice.

Storage method

- ① Place this product in its case on storage, avoiding places exposed to direct sunlight and of high temperature/humidity.

SARANCHO MODEL
CORPORATION

M175 Sakamoto Suction Simulator Instruction manual

Specifications

The following suction training techniques are possible.

- ① Suction of secretions in the oral cavity
- ② Suction of secretions in the nasal cavity
- ③ Suction of secretions in inserted tubes
- ④ Suction of secretions in the tracheotomy cannula

Contents

- ① Main body (including the bronchi, left and right lungs, and stomach)
- ② Simulated sputum
- ③ Baby powder
- ④ Funnel for rinsing
- ⑤ Tube for rinsing

Prices of replaceable parts of the Sakamoto Suction Simulator for change

Simulated sputum.....6,000 yen

※ Since this product and parts sold separately are under development, their specifications are subject to change without prior notice.

Procedure for use

- ① Before use, be sure to read <precautions>.
- ② Take the product out of the hard case and place it on a non-slip surface.
- ③ Confirm the firm connection of the bronchi/bilateral lungs/stomach.
- ④ Infuse attached, simulated sputum into the area for suction training, and initiate training.

- ⑤ After training, the bronchi can be removed and washed.

See the explanation on the accompanying sheet for washing methods and the procedure after use.

SAKAMOTO MODEL
C O R P O R A T I O N

Headquarters

34, Higashitakagicho, Shimogamo
Sakyo-ku, KYOTO 606-0865, JAPAN

Phone: +81-75-701-1135

Fax: +81-75-722-5638

E-mail: info@sakamoto-model.com

URL: <http://sakamoto-model.co.jp>

International Sales Department

Osaka factory

1-9-29, Senbahigashi, Minoh City, Osaka
562-0035, JAPAN

Phone: +81-72-729-8350

Fax: +81-72-729-8385

E-mail: trade@sakamoto-model.com