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Swallowing Mechanism Model

Learn the mechanism of swallowing, and practice practical way of preventing miss-swallowing!

This Swallowing Mechanism Model is a medical teaching equipment that was developed with the purpose of helping helpers learn the mechanism of swallowing and the proper method of preventing swallowing disorders (miss-swallowing) when giving assistance in eating. It should help prevent pneumonia caused by swallowing disorders which is often overlooked in caring for the elderly.

Practice items

- Study the mechanism of swallowing
- Learn what causes miss-swallowing
- Study proper posture when eating and its relation with bed angle
- Study neck angle (head angle) and its relation with miss-swallowing
- Practice emergency care when food gets stuck
- Swallowing drills and practice in giving mouth care

M163

Swallowing Mechanism Model

¥450,000

Supervision/ Shizuoka Junior College Social
Welfare Course Assistant Professor Prof. Ikuko Ishino

- Learn how to insert tube through nose, and the intermittent tube eating method
- Study internal mouth absorbency

Example of proper posture when eating

bed angle 30°
neck angle (head angle) 15°

Mouth care

The tongue is soft and able to pull out. Can have swallowing drills, study the tongue and muscles surrounding the mouth, and practice absorption.

Can confirm condition of the tube inserted through the nose, while observing the angle.

Can learn intermittent tube eating method, while observing the angle.

- Neck angle (head angle) can be bent about 15° forward and about 45° backward.
- Can comprehend the relation between neck angle (head angle) and movement of windpipe / Eesophagus.

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〒606-0865 京都市左京区下鴨東高木町34
TEL:075-701-1135(代) FAX:075-722-5638
E-mail:info@sakamoto-model.com
URL:<http://sakamoto-model.co.jp>